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

Copyrights: Restek Corporation

サツマイモ (Sweet Potato)

Q-sep™ QuEChERS tubes easily remove matrix interferences.

Without extract cleanup, hexadecanoic acid from the matrix obscures peaks for all the following pesticides. (example XIC chromatogram = endosulfan I)

allethrin	fenthion
buprofezin	metolachlor
cis-chlordane	myclobutanil
trans-chlordane	oxyfluorfen
chlorpyrifos	pendimethalin
cyprodinil	pentachlorothioanisole
dacthal	pirimiphos methyl
diphenamid	triadimefon
endosulfan I	triadimenol
endosulfan II	

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

Copyrights: Restek Corporation

サツマイモ (Sweet Potato)

Without cleanup, matrix masks Endosulfan I.

Peak List
1. Hexadecanoic acid
2. Endosulfan I

Column: Rxi®-5Sil MS, 20 m, 0.18 mm ID, 0.18 µm (cat.# 43602)
Sample: sweet potato spiked with pesticide mix, extracted with acetonitrile and Q-sep™ QuEChERS EN Method 15662 extraction salts, then acidified with formic acid to pH 5
Inj.: 1.0 µL splitless (hold 1 min.), 4 mm single gooseneck liner with w/wool (cat.# 22405)
Inj. temp.: 250°C
Carrier gas: helium, constant flow
Flow rate: 1.2 mL/min.
Oven temp.: 72.5°C (hold 1 min.) to 350°C @ 20°C/min.
Det: TOFMS
Transfer line temp.: 225°C
Scan range: 45-550 amu, m/z 60, 195, 197 plotted
Ionization: EI
Instrument: Agilent 6890, LECO Pegasus III



Copyrights: Restek Corporation

サツマイモ (Sweet Potato)

QuEChERS dSPE cleanup improves quantification and identification.

Peak Integration (extracted ion chromatograms)




Copyrights: Restek Corporation


サツマイモ (Sweet Potato)

Spectral Identification

QuEChERS dSPE cleanup improves mass spectral library matches.

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Table I QuEChERS extraction and cleanup procedure for pesticides from strawberries.

Sample preparation and extraction

Sample:	10g of strawberries were homogenized and placed in a 50mL PTFE centrifuge tube (cat.# 26227)
Solvent:	10mL of acetonitrile were added to homogenate Shake for 1 minute, until uniform
Salts:	4.0g MgSO ₄ (powder or granular) 1.0g NaCl 1.0g trisodium citrate dihydrate 0.5g disodium hydrogencitrate sesquihydrate] (cat.# 26213) Salts were added and vigorously shaken for 1 minute. Sample was centrifuged and the supernatant removed for cleanup. Pesticides standards (200ng/mL) were spiked in at this point.

Sample extract cleanup

QuEChERS tubes:	1mL of supernatant from the previous step was placed into several 2mL polypropylene centrifuge tubes, each containing one of the following adsorbent mixes: A. 50mg PSA + 150mg MgSO ₄ (cat.# 26124) B. 50mg PSA + 150mg MgSO ₄ + 50mg C18 (cat.# 26125) C. 50mg PSA + 150mg MgSO ₄ + 50mg GCB (cat.# 26123)
Cleanup:	Samples were shaken with the adsorbents for 30 seconds (carbon for 2 minutes), then centrifuged to produce a clear supernatant for GC/MS analysis.
Internal standard:	Pentachloronitrobenzene in a formic acid solution, pH 5.
	PSA—primary and secondary amine exchange material. GCB—graphitized carbon black

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
Copyrights: Restek Corporation

イチゴ (Strawberry)

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Table III Pesticide percent recoveries in strawberry extracts treated with C18 or GCB dSPE, relative to PSA only.

Rt (min.)	pesticide	CAS Number	action/use	classification	C18*	GCB**
9.50	dichlorvos	62-73-7	insecticide	organophosphorus	111	116
9.67	methamidophos	10265-92-6	insecticide	organophosphorus	105	107
11.75	mevinphos	7786-34-7	insecticide	organophosphorus	112	130
12.02	<i>o</i> -phenylphenol	90-43-7	fungicide	unclassified	106	97
12.14	acephate	30560-19-1	insecticide	organophosphorus	128	147
13.89	omethoate	1113-02-6	insecticide	organophosphorus	120	119
14.74	diazinon	333-41-5	insecticide	organophosphorus	108	127
14.98	dimethoate	60-51-5	insecticide	organophosphorus	124	151
15.69	chlorothalonil	1897-45-6	fungicide	organochlorine	125	13
15.86	vinclozolin	50471-44-8	fungicide	organochlorine	102	98
16.21	metalaxyl	57837-19-1	fungicide	organonitrogen	105	117
16.28	carbaryl	63-25-2	insecticide	carbamate	114	111
16.60	malathion	121-75-5	insecticide	organophosphorus	124	160
16.67	dichlofluanid	1085-98-9	fungicide	organohalogen	122	103
17.51	thiabendazole	148-79-8	fungicide	organonitrogen	88	14
17.70	captan	133-06-2	fungicide	organochlorine	88	91
17.76	folpet	133-07-3	fungicide	organochlorine	108	63

Copyrights: Restek Corporation

イチゴ (Strawberry)

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Table III Pesticide percent recoveries in strawberry extracts treated with C18 or GCB dSPE, relative to PSA only.

Rt (min.)	pesticide	CAS Number	action/use	classification	C18*	GCB**
18.23	imazalil	35554-44-0	fungicide	organonitrogen	115	95
18.39	endrin	72-20-8	insecticide	organochlorine	104	101
18.62	myclobutanil	88671-89-0	fungicide	organonitrogen	119	114
19.07	4,4-DDT	50-29-3	insecticide	organochlorine	102	95
19.22	fenhexamid	126833-17-8	fungicide	organochlorine	118	77
19.40	propargite 1	2312-35-8	acaricide	organosulfur	110	95
19.43	propargite 2	2312-35-8	acaricide	organosulfur	121	114
19.75	bifenthrin	82657-04-3	insecticide	pyrethroid	106	81
20.04	dicofol	115-32-2	acaricide	organochlorine	98	54
20.05	iprodione	36734-19-7	fungicide	organonitrogen	118	90
20.21	fenpropathrin	39515-41-8	insecticide	pyrethroid	113	96
21.32	<i>cis</i> -permethrin	52645-53-1	insecticide	pyrethroid	106	65
21.47	<i>trans</i> -permethrin	51877-74-8	insecticide	pyrethroid	109	71
23.74	deltamethrin	52918-63-5	insecticide	pyrethroid	97	52

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